

Chair Slider Belt

Description:

The padded Chair Slider Belt is a safe and comfortable way to ensure patients do not slide out of their chair, and onto the floor. Designed to encase both legs as well as the waist, the patient can move around, twisting their legs and body, without "escaping". The Chair Slider Belt is suitable for use on patients under a Restraint Authorisation.

Dimensions:

Fits Waist: Up to 130cm

Fits Thigh: Free size

The Chair Slider Belt can be hand or machine washed up to 80°C. Do up the buckles before washing and air dry.

Instructions for Use:

Ensure the patient is wearing trousers – this will protect the patient's dignity.

With the patient seated in the appropriate chair, place the belt around the patient's waist and buckle behind the chair. This rear fastening will prevent the patient from playing with the buckles. Position the leg straps under the legs, pass them through the waist loops and buckle them up behind the chair, adjusting the length as appropriate.

Check the patient regularly, particularly confirming that the straps are not placing undue pressure on the patient's thighs or pelvis. The patient should be repositioned if this is the case.

The belt is of an appropriate length to fit lounge chairs. It is recommended that "spare" length of the belt is secured away from the patient once the belt is in place.



*The Chair Slider Belt should only be used under a Restraint Authorisation.
An inappropriate or incorrectly fitted restraint can cause injury or death.*