

Adjustable Toilet Seat Raiser



Designed to assist with transfers on and off the toilet, this adjustable toilet seat raiser is easy to install on most toilets.

The seat can be adjusted to increase the height of the toilet seat by 50 mm, 100 mm or 150 mm, and includes flip up lid and lift up armrests to facilitate lateral transfers from a wheelchair or commode.

Installs over existing toilet seat hinging points, this Toilet Seat Raiser will fit most toilets.

Note: This product is not compatible with hidden or concealed hinge type toilets.



CODE	BL1097
CAPACITY	136 kg
HEIGHT	50, 100 or 150 mm
OVERALL DEPTH	480 mm
DIMENSIONS	460 mm deep x 600 mm overall width
ARMREST HEIGHT	200 mm above toilet seat
WIDTH BETWEEN ARMRESTS	480 mm
SEAT WIDTH	360 mm
WEIGHT	3.8 kg
ARTG	176831
CONFORMS TO	AS/NZS 3973:2009
WARRANTY	1 year

THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.