

# MIGHTY VOYAGER

## MIGHTY SERIES

ITEM # F-VGRPL



U.S. Patent No. D507,769S



### ABOUT

Everyone can enjoy a day at the pool, and it doesn't have to involve drilling into your deck to make it happen. Introducing the Mighty Voyager - a whole new perspective on pool accessibility. With a unique and portable design that doesn't require extensive mounting, users are able to move the lift as needed making it the perfect option for homeowners, facilities with multiple pools, and those who are looking for a convenient or less permanent solution.

This free-standing, ADA Compliant and UL Safety Certified lift boasts a 325 pound weight capacity, achieved by utilizing a sand ballast system. The 15 inch opening allows for easy filling and removal of sand when necessary, and by purchasing sand separately, a considerable amount is saved in shipping costs. Adding to the functionality of this lift is a whisper-quiet 255 degree rotation, making pool entry seamless.

The Mighty Voyager is also available in your choice of accent colors. Choose from our standard color options, or ask about customized color options.

### FEATURES: ⚙️

- ADA Compliant & UL Safety Certified
- ADA Fixing Kit included for ADA Compliance
- 325 lb (147 kg) weight capacity
- 255 degree rotation
- Flip-up arm rests
- Stainless steel construction with powder coat finish
- Limited five (5) year warranty
- Ships common carrier, some assembly required
- Made in USA

### ACCENT COLORS 🎨

Available in a variety of accent and powder coat colors



\*Blue seat standard (not a powder coat option)

### PROUDLY SOLD BY:

Affordable recreation equipment for people of all abilities.

1.888.687.3552 | [aquacreek.com](http://aquacreek.com) | [sales@aquacreek.com](mailto:sales@aquacreek.com)

# MIGHTY VOYAGER

## APPLICATIONS

This lift can work on most pool decks for in-ground and some partially raised pools and spas. US Federal ADA requirements indicate the lift needs to be mechanically affixed to the pool deck. An ADA Fixing Kit is included with the lift. To see if the Mighty Voyager works on your pool, visit our website and download the application guide, or contact our friendly sales team.

ADA Fixing Kit



For access to CADs, Revit drawings and more visit our Partner Login page on our website and become an Aqua Creek partner.

## OPTIONS / ACCESSORIES

- Rescue spine board option for safe water extractions
- Sling seat option for people with greater support needs
- Solar charging station to conserve energy
- Powder coat and accent color options to match almost any decor
- Headrest for comfortable head support
- Upper chest strap for added torso stability
- Cycle attachment to convert your lift into an aquatic bike
- Portable design for easy removal from pool side for storage
- Cover to protect your investment from the elements\*

\* Color and material options available

## PREMIUM COVER



## SOLAR CHARGER



Accessorize your lift with a cover to keep it looking like new and a solar charging station for an environmentally friendly power option.

## CUSTOM COLOR OPTIONS

Match your decor with a variety of custom color options.

