

BetterLiving

Walking Frame, Wheels and Skis

Single button release



Compact fold

A general purpose folding walking frame with front wheels and rear skis for easier handling, this lightweight aluminium walking frame can improve balance and stability when walking.

Features easily height adjustable legs with 125 mm front wheels. Single button folding mechanism for ease of transport and storage. Foam grip handles for added comfort.



CODE	BL1234
MAXIMUM USER WEIGHT	180 kg
WEIGHT	3kg
OVERALL WIDTH	640 mm
OVERALL HEIGHT	815 > 990 mm
OVERALL DEPTH	555 mm
DEPTH WHEN FOLDED	115 mm
WIDTH BETWEEN HANDLES	430 mm
WARRANTY	1 year
ARTG	176830
	Made in China

THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.



Phone 1300 734 862
support@patienthandling.com.au