



# Etac shower stools

**Edge, Easy and Smart** shower stools and accessories, designed to assist users and carers with every day routines.

**etac**<sup>®</sup>  
Creating Possibilities

# Etac Edge shower stool

The versatility of Edge makes it the perfect shower stool.

## Space saving and remarkably stable

The triangular Edge fits perfectly in a shower corner and it thereby takes up very little space. Another benefit is that the walls will provide the user support, which makes it a very safe and stable solution. Additionally the large soft ferrules ensure a stable grip on the floor. Edge is height adjustable which enables a perfect seat height for everyone.

## Pain relieving positioning

Pain relief of the hip joints can be obtained by using saddle positioning, hence having one corner facing forward when sitting. It is also possible to set the height with a slight tilt, to facilitate standing up. A soft seat is available as an accessory for extra comfort.

## Easy handling

Edge is corrosion free and designed for easy cleaning. The light weight makes it easy to move and it is stackable for easy storage. Edge comes flat packed and is simple to assemble and adjust in height without using tools.



## Accessory



## Soft Seat


A comfortable soft seat with an instantly warm feeling.


## Material:


Cellular polyethylene


Description	Colour	Height	Width	Seat width	Weight	Item no.	
Edge Shower Stool	Grey	42–57 cm 16½–22½"	52cm 20½"	45 cm 17¾"	2.6 kg 5.7 lbs	81801010	✓
Edge Shower Stool Low	Grey	34.5–42 cm 13½–16½"	52cm 20½"	45 cm 17¾"	2.4 kg 4.4 lbs	81801040	
Edge Soft Seat	Charcoal					81801132	

**Material**  
Seat: Polypropylene  
Legs: Aluminium  
Feet: TPE

 **alkali**  
pH  
min 6

 max 85°C max 185°F  
for maximum of 3 minutes

 max 130 kg max 287 lbs



Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.

# Etac Easy shower stools

Etac Easy are comfortable and safe shower stools that suits most users.

## Comfortable seating

The stools have soft contoured seats for the best comfort. The surfaces have an anti-slip pattern, which effectively drains off water so that users always sit safely.

## Increased mobility

Easy Swivel allows the user to rotate on the shower stool with little effort and with minimal strain on the body. It is convenient to clean the swivel plate since it can easily be disassembled.

## Accessory



### Soft Cushion

The soft cushion snaps easily onto the seat and is soft and warm to sit on. It's made on non slip Polypropylene



Etac Easy Swivel

Etac Easy



### Easy to fit

Practical lock function makes it easy to both fit and remove the legs.



### Easy and safe

Height adjustment has built-in safety. The push button locks against the hole when weight is put on the shower stool.



### Stability

The legs of the stool are angled outwards to increase stability. This means that the higher seat height, the greater the support base. The soft ferrules adjust perfectly to the floor.



### Comfortable and safe seats

The concave seats are comfortable to sit on. The flexibility in the leg fixture entails that the stools can handle any weight displacement which can occur when you move your upper body.

Description	Colour	Width	Depth	Seat Height	Weight	Item no.
Easy	Grey	ø 40 cm 15¾"	38cm 15"	49.5 – 53.5 cm 19½ – 21"	2.8 kg 6.2 lbs	81901010 ✓
Easy Swivel	Grey	ø 40 cm 15¾"	38cm 15"	49.5 – 53.5 cm 19½ – 21"	2.8 kg 6.2 lbs	81901050 ✓
Easy Soft Cushion	Dark Grey					80209501

<b>Construction</b> Seat: Polypropylene. Legs: Aluminium & polypropylene		max 85°C max 185°F for maximum of 3 minutes	max 150 kg max 330 lbs		Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.
--	--	---	---------------------------	--	--

# Etac Smart shower stools

Etac Smart are comfortable and safe shower stools that suits most users.

### Comfortable and stable seating

The legs are height adjustable and have soft ferrules which prevent slippage on a wet bathroom floor. The Smart shower stools can withstand weights up to 150 kg.

## Increased mobility

Smart Swivel allows the user to rotate on the shower stool with little effort and with minimal strain on the body. It is convenient to clean the swivel plate since it can easily be disassembled.



## Accessory



### Soft Cushion

Makes the seat comfortable, warm and non-slippery. Snaps on to the seat.

Material: Polypropylene



### Easy to fit

Practical lock function makes it easy to both fit and remove the legs.



## Easy and safe

Height adjustment has built-in safety. The push button locks against the hole when weight is put on the shower stool.



## Smart in the bathtub

Etac Smart can also be used in the bathtub where it provides stable and comfortable seating.

### Comfortable and safe seats

The concave seats are comfortable to sit on. The flexibility in the leg fixture entails that the stools can handle any weight displacement which can occur when you move your upper body.

Description	Colour	Width	Depth	Seat Height	Weight	Item no.
Smart	Grey	43cm 17"	38cm 15"	42 - 57cm 16½ - 22½"	3.2 kg 7.1 lbs	81951010
Smart Swivel	Grey	43cm 17"	38cm 15"	42 - 57cm 16½ - 22½"	3.2 kg 7.1 lbs	81951050
Smart Soft Cushion						80209502

## Construction

Seat: Polypropylene.  
Legs: Aluminium  
& polypropylene



max 85°C max 185°F  
for maximum of 3 minutes

max 150 kg max 330 lbs



Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.



 **etac**®  
Creating Possibilities

"Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances."

 = Standard Stock Item