

## Bed Cradle Adjustable



This height adjustable bed cradle is designed to support the weight of bedding above the legs and feet.

The lower frame slides underneath the mattress, securing the cradle to the bed, while the upper frame is positioned horizontally over the lower limbs. Bedding is then placed over the frame to remove any weight or pressure from the sensitive limbs below.

Constructed from sturdy powder-coated steel tubing. Ten different height settings to suit individual user requirements.

CODE	BL9191
HEIGHT	680 > 450 mm
WIDTH	355 mm
DEPTH	Upper 705 mm Lower 560 mm
WEIGHT	2.6 kg
WARRANTY	1 year

