

BetterLiving™

Lightweight Indoor/Outdoor Wheeled Walker

The Lightweight Indoor/Outdoor Wheeled Walker is designed for general indoor and outdoor use with a lightweight aluminium frame and 8 inch (200 mm) wheels for easy operation and manoeuvrability on smooth surfaces. The moulded handgrips with soft touch brakes are ideal for people with limited strength or mobility. Handle height can easily be adjusted, with an overall extension range of 784 - 895 mm.

Includes the added convenience of a removable shopping bag and a water resistant PVC seat for improved infection control and easier cleaning.



Specifications

NAME	Lightweight Indoor/ Outdoor Wheeled Walker
CODE	BL0060
COLOUR	Red
PRODUCT WEIGHT	7.5 kg
OVERALL HEIGHT	784 > 895mm
OVERALL WIDTH	606 mm
OVERALL DEPTH	704 mm
SEAT HEIGHT	540 mm
SEAT WIDTH	345 mm
SEAT DEPTH	31.5 mm
WHEELS	8" (200 mm)
FRAME	Aluminium
SAFE WORKING LOAD	136 kg
ARTG	176830
WARRANTY	1 Year