

## Mobile Shower Commode

This Mobile Shower Commode is designed for aged care, community and home care environments. The frame and its major components are constructed from corrosion resistant aluminium and powder coated for extra protection. All accessories, including pins, nuts and bolts are constructed with corrosion resistant materials to help ensure the strength and longevity of the product when used as intended.



### Features

- ⑥ Multi-function design for improved usability
- ⑥ Universal fitment suits most toilet types
- ⑥ Height adjustable for improved ergonomics and fitment
- ⑥ Padded accessories for user comfort
- ⑥ Corrosion resistant construction materials
- ⑥ 4x lockable castors for improved safety and stability
- ⑥ Swing back arm supports for easy transfer
- ⑥ Enclosed pan and lid with improved spill protection
- ⑥ Flip up foot rests with removable leg supports
- ⑥ Flat pack design for compact and easy storage and/or transport

### Specifications

NAME	Mobile Shower commode
CODE	BL0030
COLOUR	WHITE
PRODUCT WEIGHT	12 kg
OVERALL HEIGHT	930-1005 mm
OVERALL WIDTH	550 mm
OVERALL DEPTH	840 mm
SEAT HEIGHT	495-570 mm
SEAT WIDTH	440 mm
SEAT DEPTH	435 mm
WIDTH BETWEEN ARMRESTS	450 mm
CASTORS	4 x Lockable
FRAME	Aluminium
SAFE WORKING LOAD	136 kg
CONFORMS TO	AS/NZS 3973:2009
ARTG	176831
WARRANTY	1 Year