

BetterLiving

Shower Chair

Aluminium



Designed for users to a maximum capacity of 110 kg, this sturdy shower chair features integrated armrests for easy transfers and a comfortable moulded plastic seat with drainage slots and contoured back support.

The lightweight aluminium frame offers added durability, enhanced lifespan and increased resistance to corrosion.

The height adjustable legs include rubber tips for added grip; and each leg adjusts independently to match uneven surface requirements. One leg includes additional incremental settings to allow for more exact positioning.



THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.

CODE	BL1104
CAPACITY	110 kg
OVERALL WIDTH	540 mm
SEAT WIDTH	445 mm
BETWEEN ARMRESTS	485 mm
OVERALL DEPTH	500 mm
SEAT DEPTH	370 mm
OVERALL HEIGHT	810 > 910 mm
SEAT HEIGHT	490 > 590 mm
WEIGHT	2.9 kg
	ARTG 176831
CONFORMS TO	AS/NZS3973:2009 AS/NZS4688.2-2000 AS/NZS4688.3-2000
WARRANTY	1 year



Phone 1300 734 862
support@patienthandling.com.au