

Etac Stapel shower stool

etac[®]

Creating Possibilities

Stapel – a reliable bathroom product.



Stapel hygiene stool

Stapel is a sturdy and safe stool which can be used in many situations. It is ideal to use when stepping into the bathtub or as a footstool when using the toilet. Stapel has a tubular steel frame and a plastic top that can easily be snapped off for cleaning.



Stapel makes stepping in and out of the bathtub safe.

Description	External Dimensions	Standing Surface	Height	Item no.
Stapel	Depth 48.5 cm 19" Width 49 cm 19¼"	Depth 38 cm 15" Width 43 cm 17"	17.9 cm 7"	81204002

 alkali ph min 6	max 85°C max 185°F for maximum of 3 minutes	max 130 kg max 287 lbs		Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.
---	---	---------------------------	---	--

**iPATIENT
HANDLING™**

Phone 1300 734 862

support@patienthandling.com.au