


## Toilet Surround



Useful for transfers on and off the toilet, this portable and freestanding frame is designed to fit most toilets.

Suitable for users up to 190 kg, the sturdy steel frame is powder coated for added protection and design appeal.

The height adjustable frame features moulded anti-slip armrests for added comfort, and anti-slip rubber tips for stable grip, even on wet surfaces.



CODE	BL1165
WEIGHT	4.3 kg
CAPACITY	190 kg
OVERALL HEIGHT	680 > 810 mm
OVERALL WIDTH	580 mm
OVERALL DEPTH	440 mm
BETWEEN ARMRESTS	480 mm
	ARTG 176831
CONFORMS TO	ISO 11199-1:1999 Walking aids manipulated by both arms
WARRANTY	1 year