

# MEYRA

## AVANTI

### Economical & Optimised

#### ADAPTABILITY WITHOUT COMPROMISE

You don't have to spend a packet to get a chair optimised for you. The Avanti comes in 13 different frame colours, so make it your own. With a folding frame and removable wheels, transport is a breeze. Choose between swing-away or removable leg rest for optimal transfers!



**+** **OPTIMISED**  
DOWN TO THE  
LAST DETAIL



**Crash tested**  
Complies with international  
ISO safety standards for  
crash safety

**+** **3° wheel camber**

*Increases stability and ease of propulsion,  
without adversely affecting chair width*

**+** **Fantastic manoeuvrability**

*Lightweight and adjustable means easy  
to manoeuvre*

**+** **Pressure brake**

*Easy to activate brakes for safety and  
comfort*

**+** **Lower leg length perfectly adjustable**


*Adjust step-wise or progressively to  
perfectly fit your needs*



Adjust height-adjustable  
side guard with one hand



Telescopic and angle-  
adjustable back



Seat width in 420 mm  
and 440 mm



Compact or standard  
upper leg rest, removable  
& swing-away



Phone 1300 734 862  
support@patienthandling.com.au

# MEYRA

## AVANTI



Lower leg length can be adjusted step-wise or progressively



Up to 12° seat inclination



Schwalbe ONE tyres, hub, rim, hand rims in blue, red or yellow



Curve GEKKO, ergonomic grip, ideal for tetraplegic users

### TECHNICAL DATA

DESCRIPTION	351 SHORT FRAME	352 MEDIUM FRAME	353 LONG FRAME
Seat width	300 - 360	360 - 580	360 - 580
Overall width	Seat width + 200	Seat width + 200	Seat width + 200
Folded width	280	280	280
Seat depth	330 - 380	350 - 460	480 - 530
Seat height front	425 - 495	440 - 560	440 - 560
Seat height back	400 - 465	385 - 520	385 - 520
Overall height	830	920	920
Back height	340 - 400	380 - 500	380 - 500
Armrest height (CODE 81)	165 - 240	200 - 290	200 - 290
Length with footrests	880	1040	1120
User weight	75	135	160
Min. chair weight	12	13	14
Transport weight	8	8.5	9
Rear wheel (22")	560	560	560
Rear wheel (24")		610	610

\* with side guard CODE 100 SB + 180

Dimensions in mm, weights in kg, wheels in inches, speeds in mph (km/h), dimension tolerances  $\pm 10\text{ mm} \pm 2^\circ$ , Data subject to design changes



Phone 1300 734 862  
support@patienthandling.com.au