

# MEYRA

# NANO

## ACTIVE & UNIQUE

### RIGID FRAME AND AN OPEN HEART

The Nano is a rigid-frame wheelchair with an open, minimalist frame design. Stylish, sporty, and compact. Stable but lightweight. The compact frame fits in any car, even a convertible!

Make the Nano your own: customise frame colour, side panels, footrest and wheels.

**+** LIGHTWEIGHT  
GREAT HANDLING  
OPTIONS GALORE



**Crash tested**  
Complies with international  
ISO safety standards for  
crash safety



**+** Lightweight, great handling

*The Nano features a lightweight design with great handling*

**+** Open, rigid frame design

*Stylish, sporty, compact*

**+** Colours how you want them

*Pick the colour of the frame, side panels, footrest and wheels*

**+** Compact and easy to transport

*Lightweight and adjustable means easy to manoeuvre*



Open, rigid frame for great looks and a sporty feel



Compact - fits into your car with ease



Choose the colour of the frame, side panels, footrest and wheels



Optional aluminium castor rims and magnesium drive wheel rims

**i** PATIENT HANDLING™

Phone 1300 734 862  
support@patienthandling.com.au

# MEYRA

# NANO

## TECHNICAL DATA

	0°	3°
With wheel camber in degrees / mm	min. 490 / max. 690 SW +170	min. 540 / max. 740 SW + 240
Seat depth in mm	380 – 500	
Seat width in mm	min. 300 / max. 500	
Seat height at front in mm (FSH)	min. 450 / max. 530	
Seat height at rear in mm (RSH)	min. 350 / max. 510	
Overall height in mm	min. 600 / max. 970	
Back height in mm	min. 250 / max. 450	
Back strap height in mm	adjustable by +/- 15 each	
Frame, short, in mm	min. 760 / max. 840 approx. 170–245	
Frame, medium, in mm	min. 800 / max. 880	
Frame, long, in mm	min. 840 / max. 920	
User weight in kg	130	
Empty weight in kg	from 9.5	
Total permitted weight in kg	max. 139.5	
Transport weight in kg	<b>from 6.5</b> (without side panel, seat cushion, drive wheels)	

Dimensions in mm, weight in kg, wheels in inches, speed in km/h, size tolerance ± 10 mm, specifications subject to change



Phone 1300 734 862  
support@patienthandling.com.au