

MEYRA

XXL & XXL HD

The chair for rugged types

TOUGH AND STABLE

The Meyra XXL & XXL HD don't back down. Built to move users up to 250 kg, with a great stable and sturdy design, the XXL & XXL HD are tough. Thanks to the power tilt-in-space and recline function, as well as the choice of powered separate or centre mounted footrests, the XXL & XXL HD are built to be both tough and comfortable.



+ TOUGH AND COMFORTABLE

+ User weight up to 250 kg

The perfect bariatric wheelchair—XXL and XXL HD

+ Rear wheel drive

The classic stable drive configuration—great for indoors and outdoors

+ Built tough with lots of options

Reinforced components for strength and stability, with a wide range of standard and power options

+ LED lighting system


Safely light your way with long-lasting LEDs



Reinforced and built tough



Powerful 2 x 300W or 2 x 350W motors



Powered tilt-in-space



Long-lasting LED lighting system



Phone 1300 734 862
support@patienthandling.com.au

MEYRA

XXL & XXL HD

TECHNICAL DATA

Seat width	530 - 700	Max. obstacle height	60
Seat depth	430 - 600	Turning radius	950
Back rest height	530 - 570	Driving range with 70ah batteries	30 - 35
Leg length adjustment	370 - 510	Empty transport weight min.	approx. 130 kg
Armrest height adjustment	170 - 280	Client weight max XL / XXL HD	250 kg
Front seat height options (with /without cushion)	460 - 560 / 530 - 630	Driving speeds	6 kmh /10 kmh
Length with footrests	1,150	Motor options	2 x 350W
Length without footrests	850	Driving range max.	35 km
Width overall	670	Permitted gradient	13 %
Height without headrest	1,020 - 1,150	Electric tilt-in-space	0° to +18°
Transport length	850	Mechanical tilt-in-space	0° to +10°
Transport height minimum	650	Power backrest recline	-10° to +50°
Front wheel size	260 x 70 (10" x 2 3/4")	Mechanical backrest recline	-10° to +30°
Rear wheel size	356 x 75 (14" x 3")		

Dimensions in mm, Weights in kg, Speed in km/h, Tolerance ± 10 mm



Phone 1300 734 862
support@patienthandling.com.au