

BetterLiving™

Alternating Mattress Replacement/ Overlay



BetterLiving alternating pressure surfaces are designed to offer a cost-effective solution for the prevention and management of pressure related injuries.

Harnessing the latest technical features for improved clinical care and ease of use, BetterLiving offers both a mattress replacement and a mattress overlay option. Ideally suited to home care or aged care environments, these products can be used for people at risk to high risk of pressure injury.

ITEM	Mattress Replacement	Mattress Overlay
CODE	APMBL-R01	APMBL-L01
FEATURES	<ul style="list-style-type: none"> □ Essential features for pressure injury prevention and treatment □ Suited for patients at risk to high risk □ Replace existing mattress for superior support □ Cell in cell for enhanced user safety and comfort □ Automated care and therapy □ Easy to operate and maintain 	<ul style="list-style-type: none"> □ Feature rich overlay for pressure prevention □ Suited for patients at risk □ Used over existing mattress for basic support* □ Automated care and therapy □ Easy to operate and maintain
CAPACITY	180 kg	150 kg
CONFORMS	ARTG 289458, IEC 60601-1, IEC60601-1-2, IEC60601-1-11	
WARRANTY	2 Year	2 Year

BetterLiving by Novis

BetterLiving is a brand of cost effective yet clinically effective products designed for the Australian healthcare market. Our extensive product range includes beds and accessories, bathroom and hygiene, seating and positioning, mobility equipment and a range of aids for daily living.

As specialists in the supply of pressure care and manual transfer equipment solutions, Novis is pleased to launch our first range of affordable alternating pressure surfaces by BetterLiving.

Our focus is to offer the most efficient products at the most cost effective price – without compromise to product integrity or clinical patient care.



CPR

Handles

Cord Retention

18 Alternating Cells

Fixing Straps

Water resistant, multistretch cover

Digital Control Unit

*storage carry bag included

Why choose BetterLiving?

Safety

No more concerns over correct pressure level settings with our AutoCair feature – automatic weight sensing technology that delivers fail safe patient comfort levels irrespective of their position.

With the AutoCair feature active, the system adjusts to the appropriate pressure setting based on the patient's weight, reducing the need for carer training or monitoring. Manual adjustments can still be made for individualised therapy or comfort

Quality

While we have minimised costs to bring you the most affordable products on the market, we have never minimised on quality.

Our BetterLiving pressure surfaces conform to all the relevant Australian standards. Made using the best quality materials, every product must pass our stringent Quality Control inspection and regular testing protocols.

Simplicity

Offering all the features you need, the BetterLiving pressure surfaces are designed for basic ease of use and minimal maintenance. Easy to set up, simple to operate and straightforward to clean and maintain, these products require minimal training or support.

THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.