

## Euro Lightweight Wheeled Walker



Weighing just over 6 kilograms, the stylish Euro Lightweight Wheeled Walker is constructed from lightweight aluminium and designed for easy, every day use.

Height adjustable to best match the user, handle height is easily adjustable with an extension range of 100 mm to suit users of average height and weight. Moulded handgrips with soft touch brakes are ideal for people with limited strength.


Features a soft, wide nylon seat and stable backrest for maximum comfort when seated. Also includes a large, detachable, weatherproof fabric bag (with removable shoulder strap) that can hold up to 4.5 kg of shopping or other personal belongings; and a stick holder to ensure safe storage and transportation of walking stick.

This stylish European designer piece folds effortlessly with only one hand – simply pull the handle to collapse the walker for compact transportation or storage.



# Euro Lightweight Wheeled Walker

COLOUR	Gun Metal
CODE	BL0062
WEIGHT	6.8 kg
CAPACITY	136 kg
CONFORMS TO	ARTG 176830 and ISO 11199-2:2005
WARRANTY	1 Year



THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.