

Air Transfer Mat

Float on a cushion of air.

Safer and easier lateral transfer for patients and carers.




**REINFORCED
HANDLES**



**DUAL CONNECTION
POINTS**



**COMPATIBLE
PUMP DESIGN***



**SWL UP
TO 450KG**

Air Transfer Mat

GREATER STAFF SAFETY AND PATIENT COMFORT

- 80-90% less physical effort needed to perform transfers, greatly reducing injury risks to carers.
- The Novis Air Transfer Mat cradles and supports the patient, maintaining a stable position while transferring.
- Proven design eliminates shearing and friction forces providing more comfortable patient transfers.

AVAILABLE IN 3 SIZES

CODE	SIZE (W x L)
ATM102D	864x1980 mm
ATM103D	990x1980 mm
ATM104D	1270x1980 mm

BLOWER PUMP

ATMP001	Air Transfer Mat Blower Pump
---------	------------------------------

*All Novis Air Transfer Mats are compatible with many of the most commonly available pump units

INFECTION CONTROL AND COST EFFECTIVENESS

- Single patient use mats deliver greater infection control management.
- Save time and costs on decontamination as no laundering is required.
- Greatly reduce costs associated with reductions in caregiver/nursing injuries.
- Number of staff required is significantly reduced during manual handling transfers representing cost and time efficiency.

EFFECTIVE FOR A WIDE RANGE OF PATIENT TRANSFER & MOVEMENT REQUIREMENTS

- Accommodates various patient sizes up to 450kg.
- Vapour permeable nylon helps manage microclimate while helping provide added protection against shear and friction.
- Heavy duty patient safety straps and handles.
- Radiolucent/MRI safe.



1. Electrical Power Cord
2. Support Brackets
3. Pump Controls On/Off
4. Hose Connection Point
5. 2nd Hose Connection Point
6. Colour Coded Transfer Handles
7. Safety Strap Buckles (only fully tighten when inflated)
8. Air Hose
9. Cord Retention Strap