

CairMax Duo

A versatile foam air support surface with the ability to alternate between static or dynamic support, in response to the needs of the client. With a focus on support and versatility, the mattress functions as a static pressure relieving surface without the control unit, however when clinically required, can “step up” to provide dynamic support with the control unit connected.



Versatility

Mattress and control unit can be purchased as a system or separately – place a CairMax Duo mattress on every bed frame to provide all clients with superior static support then add a control unit for increased therapy as needed.

Simplicity

This ‘plug and play’ system is very simple to use and requires minimum training to operate, making it ideal for community or aged care use.

Safety

No more concerns over correct pressure level settings with our AutoCair feature – automatic weight sensing technology that automatically adjusts the pressure settings to keep the client within a safe pressure range for effective therapy.

Comfort

A durable foam core combined with alternating air cells embedded with castellated foam, and a memory foam topper, delivers comfort and support for clients up to 230 kilograms.



CairMax Duo	Standard	2000 x 880 x 180 mm	FAMCM-R01
System	King Single	2000 x 1050 x 180 mm	FAMCM-R01K
CairMax Duo,	Standard	2000 x 880 x 180 mm	FAMCM-M01
Mattress Only	King Single	2000 x 1050 x 180 mm	FAMCM-M01K
CairMax Duo, Control Unit Only			APMBL-CU01



Phone 1300 734 862
support@patienthandling.com.au