

# scian Wrist Blood Pressure Monitor

Simultaneously displays systolic pressure,  
diastolic pressure and pulse rate

Lightweight, portable design

Automatically inflates to correct pressure

Single button operation



This convenient, portable blood pressure monitor wraps around the wrist and uses the oscillometric principle to simply and accurately measure systolic pressure, diastolic blood pressure and pulse rate.

With an easy fitting wrist cuff and one touch operation, the unit automatically inflates to the correct pressure, using oscillometric technology to simultaneously measure the arterial blood pressure and pulse rate.

A 90 memory storage capacity allows the review of the last 90 readings.

Also features low battery indicator and automatic power down after 3 minutes of no operation. Requires 2 x AAA batteries (included) and suits wrist circumference of 125 to 205 mm.

The compact unit is lightweight and easy to carry when travelling or for daily use.

<b>CODE</b>	BSSDS733
<b>CUFF</b>	125 > 205 mm circumference
<b>HEIGHT</b>	25 mm
<b>WIDTH</b>	61 mm
<b>DEPTH</b>	62 mm
<b>WEIGHT</b>	190 g with batteries
<b>POWER</b>	2 x AAA 1.5V batteries
<b>ACCURACY</b>	Blood Pressure $\pm 3$ mmHg Pulse Rate $\pm 5\%$
<b>CONFORMS TO</b>	ARTG 165748
<b>WARRANTY</b>	1 year

FOR PEOPLE WITH HIGH BLOOD PRESSURE. CONSULT YOUR MEDICAL PROFESSIONAL TO EVALUATE THE READINGS. CHECK YOUR DEVICE PERIODICALLY FOR ACCURACY.



Phone 1300 734 862  
support@patienthandling.com.au