

Lowback Backrest Support



Fabric: open weave, non-allergenic brushed nylon


Core: three grades of high-resilient urethane foam

Frame: impact resistant polycarbonate biomechanically moulded



Transforms ordinary chairs into ergonomically correct seating

Removable, adjustable lumbar pad



An award-winning portable backrest that transforms ordinary chairs into ergonomically correct seating to provide relief from back pain.

Unlike ordinary "L" shaped chairs, the "S" shape of this backrest moulds the spine into a more anatomically correct position, helping to enhance overall posture by supporting proper spinal alignment.

- lightweight and portable; includes elastic positioning straps and handle
- unique polycarbonate frame with polyurethane foam gently support natural spine contours

- removable and adjustable lumbar pad for additional, customisable back support
- hypoallergenic brushed nylon cover can be removed and sponge washed
- helps prevent fatigue, provide relief from pain aggravated by poor posture and ease pain caused by spinal disc problems and sciatica

CODE	OBUSLB
COLOUR	Black
WIDTH	450 mm
HEIGHT	530 mm
DEPTH	75 mm
MATERIALS	Nylon, Foam, Polycarbonate
WEIGHT	1kg
CAPACITY	Unlimited
WARRANTY	1 Year

ENDORSED BY



Canadian Physiotherapy Association

