

3501 & 3508 Pagers

cura¹
Helping To Care



3508

Monitor up to three
cordless devices



3501

Monitor up to nine
cordless devices




***... easy cordless monitoring
for facility & home care***

3501 & 3508 Pagers

Wireless Pagers are ideal for carers on the move and their wide compatibility with the cura¹ cordless range makes them impressively versatile. The pagers boast advanced features and reliability yet are simple and intuitive to use.


3501 ACTIVECARE LCD PAGER

- Large easy-to-read 7-digit LCD screen displays device type and ID.
- Monitor up to 9 wireless devices.
- Up to 5 simultaneous alarms can be stored in the pager's "Call Review" memory.
- Vibration and tone selection for hard of hearing.
- Out of range check.




3508 LED PAGER

- Can monitor up to 3 cordless devices.
- Vibration and tone selection for hard of hearing.
- Auto-reset when paired with a 2719 Transmitter.
- Direct signalling from cura¹ cordless devices for use in home care environments.
- Easily cleanable buttons for a hygienic setting.



WHAT'S IN THE BOX?



3501 or 3508 Pager

Batteries

QuickStart User Guide

COMPATIBLE PRODUCTS

- All cura¹ cordless pads, mats, buttons & sensors.
- 3506 Soft Cover
- 3504 Mounting Bracket
- 3010 SeizureAlert Monitor



Certifications may vary between products.

cura¹
Helping To Care

10 Gray Street, Kilkenny, South Australia 5009

sales@cura1.com | www.cura1.com

+61 (0) 8 8268 8343 (International)

1300 12-CARE (2273)

PB-003