

The Compact

Easy to lift, fold
and go



Lightweight Compact Wheeled Walker

Weighing just 7.2 kilograms, the stylish Compact Wheeled Walker is constructed from lightweight aluminium and designed for easy, every day use.

Height adjustable to best match the user, handle height is easily adjustable with an extension range of 120 mm to suit users of average height and weight. Moulded handgrips with soft touch brakes are ideal for people with limited strength.

Features a soft, wide nylon seat and stable backrest for maximum comfort when seated. Also includes a large, detachable, weatherproof fabric bag (with removable shoulder strap) that can hold up to 5 kg of shopping or other personal belongings; and a stick holder to ensure safe storage and transportation of walking stick.

This stylish designer piece folds effortlessly with only one hand – simply pull the handle to collapse the walker, secure the seat clasp and then fold the handles down towards the wheels for double folding. The wheeled walker easily folds and can be carried without any extra securing needed, making it ideal for transport and storage.



The Compact


Easy to lift, fold and go



COLOUR		Matte Black
CODE		BL00012
WEIGHT	FRAME	7.2 kg
	BAG AND BACKREST	Additional 0.4 kg
CAPACITY		136 kg
CONFORMS TO		ISO 11199-2:2021
WARRANTY		1 Year



HEIGHT




820 mm min
940 mm max

LENGTH




696 mm

WIDTH




620 mm

SEAT DEPTH




250 mm

SEAT WIDTH




458 mm

SEAT HEIGHT



545 mm

FOLDED WIDTH



230 mm

CASTOR SIZE



200 mm

THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.