

# The Premium

Comfort and class  
in every step



## Premium Lightweight Wheeled Walker

Weighing just 6.5 kilograms, the stylish Lightweight Wheeled Walker is constructed from lightweight aluminium and designed for easy, every day use.

Height adjustable to best match the user, handle height is easily adjustable with an extension range of 120 mm to suit users of average height and weight. Moulded handgrips with soft touch brakes are ideal for people with limited strength.

Features a soft, wide nylon seat and stable backrest for maximum comfort when seated. Also includes a large, detachable, weather resistant fabric bag (with removable shoulder strap) that can hold up to 5 kg of shopping or other personal belongings; and a stick holder to ensure safe storage and transportation of walking stick.

This stylish designer piece folds effortlessly with only one hand – simply pull the handle to collapse the walker for compact transportation or storage.



# The Premium


Comfort and class  
in every step



COLOUR		Green
CODE		BL00011
WEIGHT	FRAME	6.5 kg
	BAG AND BACKREST	Additional 0.7 kg
CAPACITY		136 kg
CONFORMS TO		ISO 11199-2:2021
WARRANTY		1 Year



HEIGHT




850 mm min  
970 mm max

LENGTH




705 mm

WIDTH




610 mm

SEAT  
DEPTH




250 mm

SEAT  
WIDTH



461 mm

SEAT  
HEIGHT



545 mm

FOLDED  
WIDTH



230 mm

CASTOR  
SIZE



250 mm  
& 200 mm

THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.